



Shelf 1: Height 4' 5.00", Length 6' 0.00", Depth 1' 8.00"

Large Planogram

UPC	Johnstone Part #	Fluke Model #	Description	Min Units	Suggested List
9596979446	H28-271	FLUKE-279FC/IFLEX	TRUE RMS THERMAL MULTIMETER WITH IFLEX	1	\$1,199.99
9596972419	H85-825	FLK-3000FC HVAC	FC WIRELESS HVAC SYSTEM	1	\$629.99
9596906180	H25-862	FLUKE-C0-220	CARBON MONOXIDE METER	1	\$459.99
9596924524	H26-367	FLUKE-87-5	INDUSTRIAL TRUE RMS MULTIMETER WITH TEMPERATURE	1	\$429.99
9596907133	H25-943	FLUKE-52-2	DUAL INPUT THERMOMETER	1	\$355.99
9596928751	H26-513	FLUKE-971	TEMPERATURE HUMIDITY METER	2	\$349.99
9596982693	H28-396	FLUKE-180LR	2 LINE LASER LEVEL	2	\$339.99
9596979834	H28-203	FLUKE-902FC	TRMS WIRELESS HVAC CLAMP METER	2	\$299.99
9596965920	H27-621	FLUKE-116/323 KIT	HVAC MULTIMETER AND CLAMP METER COMBO KIT	1	\$279.99
9596985404	H28-521	T6-1000	ELECTRICAL TESTER WITH FIELDSENSE, 1000V	2	\$249.99
9596933253	H26-578	FLUKE-561	INFRARED THERMOMETER FOR THE HVAC PROFESSIONAL	2	\$202.99
9596962375	H27-469	FLUKE-324	400A AC TRUE RMS CLAMP METER WITH TEMPERATURE	2	\$199.99
9596906181	H25-864	FLUKE-CO-205	ASPIRATOR KIT	2	\$199.99
9596904077	H25-585	80PK-8	TYPE-K, THERMOCOUPLE PIPE CLAMP PROBE ACCESSORY	2	\$181.99
9596932419	H25-748	FLUKE-116	HVAC TRUE RMS MULTIMETER WITH TEMPERATURE	2	\$179.99
9596985411	H28-520	T6-600	ELECTRICAL TESTER WITH FIELDSENSE, 600V	2	\$179.99
9596961355	H27-464	FLUKE-414D	LASER DISTANCE METER 50M/165FT	2	\$142.99
9596962020	H27-462	FLUKE-62 MAX	10:1 NON-CONTACT IR THERMOMETER	3	\$101.99
9596969071	H27-700	FLUKE-59 MAX NA	8:1 NON-CONTACT INFRARED THERMOMETER	3	\$81.99
9596920645	H26-483	TL220	SUREGRIP INDUSTRIAL TEST LEAD SET	2	\$81.99
9596948492	H26-925	RLD2	UV REFRIGERANT LEAK DETECTOR AND FLASHLIGHT	2	\$59.99
9596936436	H26-654	80PK-11	TYPE-K, VELCRO TEMPERATURE PROBE	2	\$45.99
9596907766	H25-937	TPAK	METER HANGING KIT	2	\$45.99
9596937096	H26-059	80BK-A	TYPE-K, INTEGRATED DMM PROBE	2	\$45.99
9596900637	H24-758	C25	SOFT CARRYING CASE	2	\$41.99
9596900974	H24-751	80PK-1	TYPE-K, THERMOCOUPLE MODULE W/BEAD PROBE	3	\$35.99
9596905719	H25-839	H5	HOLSTER FOR THE T SERIES TESTERS	2	\$35.99
9596952154	H27-196	TL175	TWISTGUARD TEST LEADS	3	\$34.99
9596937097	H25-749	80AK-A	TYPE-K, THERMOCOUPLE ADAPTER	3	\$30.99
9596961242	H27-554	AC175	THREADED ALLIGATOR CLIP SET	3	\$30.99
9596929078	H25-736	FLUKE-1AC-A1-II	AC NON-CONTACT VOLTAGE DETECTOR, 90V - 1000V	12	\$29.99
9596929087	H25-863	FLUKE-1LAC-A-II	AC NON-CONTACT LOW VOLTAGE DETECTOR, 20V - 90V	12	\$29.99
9596956442	H27-195	FLK2A/90-1000VC	AC NON-CONTACT VOLTAGE DETECTOR, 90-1000V	12	\$29.99
9596901158	H24-756	TL75	HARD POINT TEST LEAD SET	3	\$25.99
9596905769	H25-840	TP1	TEST PROBE FLAT BLADE	2	\$20.99
9596904599	H25-747	203403	FUSE, 11A/1000V LOT QTY 1	2	\$14.40
9596904601	H25-746	203411	FUSE, 440mA/1000V LOT QTY 1	2	\$9.60