

Meters & Testers

Fluke

FLU-TB-B-1

Originated: 02-20-04 Revised: 9-11-17

EXAMPLE:
Peg: (1, 20-21)
References the row
and column
designation for the
peghooks. Row: 1,
Column: 20-21

Revision: Plan-o-gram includes new
Fluke T6 Electrical Testers with
FieldSense.

Use this space to display Fluke products not on this planogram.		Wire Lit Holder : (57, 2-8) Fluke Catalog		4" Peg: (4, 3-4) TL223		4" Peg: (3, 3-4) FLUKE-180LR		4" Peg: (2, 3-4) FLUKE-376 FC		4" Peg: (1, 4-5) FLUKE-87-5	
				4" Peg: (4, 10-11) TL220		4" Peg: (3, 9-10) TPAK		4" Peg: (2, 9-10) FLUKE-374 FC		4" Peg: (1, 10-11) FLK-3000 FC	
		4" Peg: (5, 11-12) TL75		4" Peg: (4, 15-16) TL175		4" Peg: (3, 15-16) T+PRO		4" Peg: (2, 15-16) FLUKE-324		4" Peg: (1, 17-18) FL UKE-117	
				4" Peg: (5, 21-22) CAMO-C25		4": (3, 21-22) LVD2.CT		4" Peg: (2, 21-22) T6-1000		4" Peg: (1, 24-25) FLUKE-115	
		4" Peg: (5, 27-28) C25				4" (4, 23-24) TL224		4": (3, 25-26) FLUKE-1AC-A1-II		4" Peg: (2, 28-29) T6-600	
				4" Peg: (5, 33-34) H5		4" : (4, 28-29) AC220		4": (3, 29-30) FLUKE-1AC-A1-II		4" Peg: (1, 31-32) FLUKE-62MAX	
		4" (4, 33-34) AC285				4": (3, 33-34) FLK2AC/90-1000VVC		4" Peg: (2, 33-34) T5-600 U4A			



LOZIER PARTS LIST

Qty	Part #	Description
30	HP1314	4" Peghook
01		Header Sign & Literature Holder (Ask your Fluke Rep.)

Revision
Literature Order # 1631613


Graybar
works to your advantage.

Planogram Item List

PLANO CODE	PFX	PURCH GRP	PART.	MATERIAL	DESCRIPTION	LOC CODE	ORD QTY
FLU-TB-B-1	FLU	PA	AC220	22052303	SUREGRIP ALLIGATOR CLIPS	A	2
FLU-TB-B-1	FLU	PA	AC285	22052855	SUREGRIP LARGE JAW ALLIGATOR CLIPS	A	2
FLU-TB-B-1	FLU	PA	AC175	25475011	THREADED ALLIGATOR CLIP SET 1000V	A	2
FLU-TB-B-1	FLU	PA	CAM0-C25	25442354	CAMOUFLAGE CARRYING CASE	A	2
FLU-TB-B-1	FLU	PA	C25	90006432	SOFT CASE	A	2
FLU-TB-B-1	FLU	PA	FLUKE-117	25079181	ELECTRICIANS TRUE RMS MULTIMETER	A	2
FLU-TB-B-1	FLU	PA	FLUKE-115	25079187	TRUE RMS DIGITAL MULTIMETER	A	2
FLU-TB-B-1	FLU	PA	FLK-3000 FC	25643667	FC WIRELESS DIGITAL MULTIMETER	A	1
FLU-TB-B-1	FLU	PA	FLUKE-324	25447164	MID SIZE TRUE RMS CLAMP METER WITH TEMPERATURE	A	2
FLU-TB-B-1	FLU	PA	FLUKE-374FC	25762587	600A AC/DC TRMS WIRELESS DIGITAL CLAMP METER	A	1
FLU-TB-B-1	FLU	PA	FLUKE-376FC	25762585	1000A AC/DC TRMS WIRELESS DIGITAL CLAMP METER W/ IFLEX PROB	A	2
FLU-TB-B-1	FLU	PA	FLUKE-180LR	25850945	2 LINE LASER LEVEL RED	A	2
FLU-TB-B-1	FLU	PA	FLUKE-62 MAX	25437751	10:1 NON-CONTACT IR THERMOMETER	A	3
FLU-TB-B-1	FLU	PA	FLUKE-87-5	22086718	IND TRUE RMS MULTIMETER TEMP	A	1
FLU-TB-B-1	FLU	PA	H5	97204656	HOLSTER BELT ACCESSORY	A	2
FLU-TB-B-1	FLU	PA	FLK2AC/90-1000VC	25361862	ACV DETECTOR 90- 1000V VOLTALERT COUNTER TOP	A	6
FLU-TB-B-1	FLU	PA	FLUKE-1AC-A1-II	96051655	ACV DETECTOR 90-1000V	A	6
FLU-TB-B-1	FLU	PA	LVD2/CT	25091146	VOLT DETECTOR W/LIGHT	A	6
FLU-TB-B-1	FLU	PA	T+PRO	25092665	ELECTRICAL TESTER W/GFCI TRIP	A	2
FLU-TB-B-1	FLU	PA	T6-1000	25964018	ELECTRICAL TESTER WITH FIELDSENSE, 1000V	A	2
FLU-TB-B-1	FLU	PA	T6-600	25964020	ELECTRICAL TESTER WITH FIELDSENSE, 600V	A	2
FLU-TB-B-1	FLU	PA	T5-600-USA	97204655	ELECTRICAL TESTER	A	2
FLU-TB-B-1	FLU	PA	TL175	25373911	TWISTGUARD TEST LEAD SET	A	2
FLU-TB-B-1	FLU	PA	TL220	22052860	SUREGRIP IND TEST LEAD SET	A	2
FLU-TB-B-1	FLU	PA	TL223	22052305	ELECTRICAL TEST LEAD SET	A	2
FLU-TB-B-1	FLU	PA	TL224	22052304	SILICONE TEST LEAD	A	2
FLU-TB-B-1	FLU	PA	TL75	91011135	TEST LEAD SET	A	2
FLU-TB-B-1	FLU	PA	TP175	25674927	TWISTGUARD TEST PROBES	A	2
FLU-TB-B-1	FLU	PA	TPAK	99502059	METER HANGING KIT	A	2
			NEW PRODUCTS				