



Shelf 1: Height 5' 5.00", Length 4' 0.00", Depth 1' 8.00"

UPC	Fluke Model #	Description	Min Units	Suggested List Price
9596979446	FLUKE-279FC/IFLEX	TRUE-RMS THERMAL MULTIMETER WITH IFLEX	1	\$1,199.99
9596979662	FLUKE-1587FC	2-1 INSULATION MULTIMETER	1	\$799.99
9596982692	FLUKE-180LG	2 LINE LASER LEVEL GREEN	1	\$479.99
9596979839	FLUKE-376FC	1000A AC/DC TRMS WIRELESS CLAMP METER WITH IFLEX	2	\$449.99
9596924524	FLUKE-87-5	INDUSTRIAL TRUE RMS MULTIMETER WITH TEMPERATURE	1	\$429.99
9596979846	FLUKE-375FC	600A AC/DC TRMS WIRELESS CLAMP METER	1	\$399.99
9596982693	FLUKE-180LR	2 LINE LASER LEVEL RED	1	\$339.99
9596985404	T6-1000	ELECTRICAL TESTER WITH FIELDSENSE, 1000V	3	\$249.99
9596982688	FLUKE-3PR	3 POINT LASER	1	\$239.99
9596962375	FLUKE-324	400A AC TRUE RMS CLAMP METER WITH TEMPERATURE	2	\$199.99
9596932420	FLUKE-117	ELECTRICIAN'S TRUE RMS MULTIMETER	2	\$199.99
9596984009	FLUKE-64 MAX	MULTIFUNCTIONAL IR THERMOMETER 20:1 DS	3	\$199.99
9596985411	T6-600	ELECTRICAL TESTER WITH FIELDSENSE, 600V	3	\$179.99
9596905311	T5-1000 USA	ELECTRICAL TESTER, VOLTAGE, CONTINUITY AND CURRENT, 1000V	2	\$159.99
9596977298	PRV240	PROVING UNIT PORTABLE BATTERY POWERED VOLTAGE SOURCE	2	\$149.99
9596961355	FLUKE-414D	LASER DISTANCE METER 50M/165FT	2	\$142.99
9596905312	T5-600 USA	ELECTRICAL TESTER, VOLTAGE, CONTINUITY AND CURRENT, 600V	2	\$119.99
9596932792	T+PRO	ELECTRICAL TESTER WITH GFCI TEST	3	\$99.99
9596962020	FLUKE-62 MAX	10:1 IR THERMOMETER	3	\$99.99
9596920645	TL220	SUREGRIP INDUSTRIAL TEST LEAD SET	2	\$81.99
9596910339	AC220	SUREGRIP ALLIGATOR CLIPS	2	\$45.99
9596908263	TPAK	METER HANGING KIT WITH MAGNET	2	\$45.99
9596900637	C25	SOFT CARRYING CASE	2	\$40.99
9596951635	LVD2/CT	NON-CONTACT AC VOLTAGE DETECTOR AND LED FLASHLIGHT	12	\$39.99
9596940008	C115	SOFT CASE	2	\$35.99
9596905719	H5	HOLSTER FOR T-SERIES ELECTRICAL TESTERS	2	\$35.99
9596910344	TL224	SUREGRIP SILICONE TEST LEADS	2	\$30.99
9596952154	TL175	TWISTGUARD TEST LEADS	3	\$30.99
9596961242	AC175	THREADED ALLIGATOR CLIP SET	3	\$30.99
9596956442	FLK2AC/90-1000VC	AC VOLTAGE DETECTOR 90 - 1000VAC	12	\$29.99
9596929078	FLUKE-1AC-A1-II	AC VOLTAGE DETECTOR 90 - 1000VAC	12	\$29.99
9596952155	TP175	TWISTGUARD TEST PROBES	3	\$25.99
9596905769	TP1	SLIM REACH FLAT BLADE TEST PROBES	2	\$20.49