



Fluke Corporation**Suggested Product Mix****September 11, 2017**

UPC	Fluke Model #	Description	Min Units	Suggested List Price
9596979839	FLUKE-376FC	1000A AC/DC TRMS WIRELESS CLAMP METER	2	\$449.99
9596924524	FLUKE-87-5	INDUSTRIAL TRUE RMS MULTIMETER WITH TEMPERATURE	1	\$429.99
9596985404	T6-1000	ELECTRICAL TESTER WITH FIELDSENSE, 1000V	2	\$249.99
9596962375	FLUKE-324	400A AC TRUE RMS CLAMP METER WITH TEMPERATURE	2	\$199.99
9596932420	FLUKE-117	ELECTRICIAN'S TRUE RMS MULTIMETER	2	\$199.99
9596985411	T6-600	ELECTRICAL TESTER WITH FIELDSENSE, 600V	2	\$179.99
9596905312	T5-600 USA	ELECTRICAL TESTER, VOLTAGE, CONTINUITY AND CURRENT, 600V	2	\$119.99
9596962020	FLUKE-62 MAX	10:1 NON-CONTACT INFRARED THERMOMETER	3	\$99.99
9596932793	T+PRO	ELECTRICAL TESTER WITH GFCI	3	\$99.99
9596920645	TL220	SUREGRIP INDUSTRIALTEST LEAD SET	2	\$81.99
9596908263	TPAK	METER HANGING KIT WITH MAGNET	2	\$45.99
9596900637	C25	SOFT CARRYING CASE	2	\$40.99
9596951635	LVD2/CT	NON-CONTACT AC VOLTAGE DETECTOR AND LED FLASHLIGHT	12	\$39.99
9596905719	H5	HOLSTER FOR T-SERIES ELECTRICAL TESTERS	2	\$35.99
9596952154	TL175	TWISTGUARD TEST LEADS	3	\$30.99
9596956442	FLK2AC/90-1000VC	AC VOLTAGE DETECTOR 90 - 1000V	12	\$29.99
9596929078	FLUKE-1AC-A1-II	AC VOLTAGE DETECTOR 90 - 1000V	12	\$29.99
9596901158	TL75	HARD POINT TEST LEAD SET	3	\$25.99
9596913320	TP220	SUREGRIP TEST PROBES	2	\$20.49
9596905769	TP1	SLIM REACH FLAT BLADE TEST PROBES	2	\$20.49