

SEE IT. SAVE IT. SHARE IT.

FLUKE CONNECT

FLUKE

HVAC Multimeter with Temperature

118

902 FC

FLUKE CONNECT

True-rms Clamp Meter

324

Electrical Tester

TS-400

Temperature Humidity Meter

871

Non-contact Thermometer

42 5000

TYPE K THERMOCOUPLE PPE CLAMP TEMPERATURE PROBE

BCPK-8

MultiMeter

100

Camouflage Carrying Case

5000-000

Test & Communication Adapter

3000-000

Motor Hanging Kit with Magnet and Straps

TRAK

TwistGuard Test Leads

FLUKE

Right Backprobe Pin Set

TP08

GENERAL PURPOSE BIAS PROBE

BCPK-1

Shelf 1: Height 4' 5.00", Length 2' 0.00", Depth 1' 8.00"

Small HVAC/R Plan-o-gram

Fluke Corporation**Suggested Product Mix****February 1, 2018**

UPC	Trane P/N	Fluke Model #	Description	Min Units
9596928751	FLUKE-971	FLUKE-971	TEMPERATURE AND HUMIDITY METER	1
9596979834	MTR01585	FLUKE-902FC	TRMS WIRELESS HVAC CLAMP METER	2
9596962375	FLUKE-324	FLUKE-324	400A AC TRUE RMS CLAMP METER WITH TEMPERATURE	2
9596904077	80PK-8	80PK-8	TYPE-K, THERMOCOUPLE PIPE CLAMP PROBE ACCESSORY	2
9596932419	FLUKE-116	FLUKE-116	HVAC TRUE RMS MULTIMETER WITH TEMPERATURE	2
9596985411	MTR01668	T6-600	ELECTRICAL TESTER WITH FIELDSENSE, 600V	2
9596962020	FLUKE-62 MAX	FLUKE-62 MAX	10:1 NON-CONTACT IR THERMOMETER	3
9596985447	CAS02144	CAMO-C25/FO	CAMOUFLAGE CARRYING CASE REAL TREE	2
9596907766	TPAK	TPAK	METER HANGING KIT	2
9596900974	80PK-1	80PK-1	TYPE-K, THERMOCOUPLE MODULE W/BEAD PROBE	3
9596952154	TL175	TL175	TWISTGUARD TEST LEADS	3
9596937097	80AK-A	80AK-A	TYPE-K, THERMOCOUPLE ADAPTER	3
9596929078	FLUKE-1AC-A1-II	FLUKE-1AC-A1-II	AC NON-CONTACT VOLTAGE DETECTOR 90 - 1000VAC	12
9596926128	TP88	TP88	RIGID BACKPROBE PIN SET	2