



Datapaq Insight Food Tracker HACCP Report	File : Batch Oven.paq	Company: Datapaq Limited	
Site : Cambridge UK	Process : Batch Oven	Product: Duck	

Probe	Fo Results		Chill		Maximum / Minimum							Time at Temperature			
			chill 1 (60.0 To 30.0°C)	Chill 2 (26.7 To 4.4°C)											
	Fo (min)	Decimal Reductions	Fall Time (mm:ss.tt)	Fall Time (mm:ss.tt)	Maximum (°C)	Max. Reached (mm:ss.tt)	Mean (°C)	Deviation From 63.0°C	Standard Deviation	Minimum (°C)	Min. Reached (mm:ss.tt)	Time Above 74.0°C (mm:ss.tt)	Time To Reach 74.0°C (mm:ss.tt)	Time Above 80.0°C (mm:ss.tt)	Time To Reach 80.0°C (mm:ss.tt)
#1 (°C) Rack 8 Left Front Edge - Air	***	***	08:35.00	20:00.00	169.1	126:40.00	43.6	+106.1	65.8	-1.0	249:15.00	120:35.00	06:30.00	120:25.00	06:35.00
#2 (°C) Rack 4 Middle Front Edge - Product	46530.67	9269.06	32:40.00	110:20.00	88.4	129:15.00	28.0	+25.4	30.4	-0.2	04:05.00	64:00.00	83:55.00	39:10.00	103:40.00
#3 (°C) Rack 6 Right Centre - Air	***	***	08:25.00	09:20.00	172.4	126:10.00	45.1	+109.4	69.0	-1.5	363:05.00	120:50.00	06:25.00	120:45.00	06:25.00
#4 (°C) Rack 6 Left Centre - Product	151688.46	30216.83	29:55.00	75:25.00	92.6	128:15.00	28.0	+29.6	32.0	-0.2	02:20.00	70:45.00	74:30.00	51:05.00	90:10.00
#5 (°C) Rack 1 Left Front Edge - Product	314575.61	62664.46	38:25.00	101:30.00	93.9	128:15.00	31.5	+30.9	33.1	-0.2	03:25.00	86:30.00	64:25.00	68:40.00	77:25.00
#6 (°C) Rack 1 Right Front Edge - Product	308310.21	61416.38	34:05.00	113:35.00	94.1	127:40.00	29.3	+31.1	32.3	-0.1	02:40.00	74:50.00	69:35.00	59:55.00	80:30.00

Probe	Slopes	Peak Difference			Area Under Curve
	Positive Slope (°C/min)		Peak Difference (°C)	Time Reached (mm:ss.tt)	Area (°C)min
#1 (°C) Rack 8 Left Front Edge - Air	220.20	 2.6 143.3 140.7 13:55.00	140.7	13:55.00	19868.42
#2 (°C) Rack 4 Middle Front Edge - Product	12.60				12986.34
#3 (°C) Rack 6 Right Centre - Air	175.20				20539.49
#4 (°C) Rack 6 Left Centre - Product	19.20				13016.68
#5 (°C) Rack 1 Left Front Edge - Product	27.60				14543.65
#6 (°C) Rack 1 Right Front Edge - Product	55.20				13598.75